QUICK START GUIDE

Non-contact Infrared Thermometer

1. Product Structure

1.Probe

2.LCD Display

3.Memory Button/Adjustable Button

4.Memory Button/Adjustable Button

5.Mode Button (Body/Surface)

6.Measure Button

7.Battery Cover

2. Display Description

1.Surface temperature

2.Body temperature

3.Display Value

4.Storage Location

5.Battery Symbol

6.Buzz Symbol

7.Temperature Unit

8.Storage Data Readout

3. How to Set Temperature Unit & Mute/Unmute Mode?

Press the “Mode” button to enter the F1 mode. Under the F1 mode you can press the +/- button to switch between °F/°C. Then continue to press the set button 4 times to exit and save the setting.

Press the “Mode” button to enter the F4 mode. Under the F4 mode, you can press the +/- button to turn on/off the beep sound. Then continue to press the set button 1 time to exit and save the setting. In the mute mode, the 🔇 will not display on the screen.

Under F1 Function

Press the +/- button to switch between °F/°C

Under F4 Function

Press the +/- button to turn on/off the beep sound

4. How To Read History Readings?

- 1.Keep the thermometer turned on.**
- 2.Press Memory Button to get into the memory mode.**
- 3.Press the +/- button again to check the 32 sets of memories one by one.**

Press Memory Button

Page 2

5. Notice

- 1.Continuous measurement will heat up the sensor, which will make the readings inaccurate. To get the accurate temperature, please wait about 10 seconds between each measurement.**
- 2.Keep a distance of 3–5 cm (1.2–2 inches) when measure forehead temperature.**
- 3.Take the forehead temperature in the correct place: The center of forehead.**

Far away from forehead

Distance of 3–5 cm (1.2–2 inches)

6. Measure Forehead Temperature

1.Press the measure button to turn on the machine.

Power on

2.Press the mode button to choose body mode.

Choose mode

3.Point the probe at the center of forehead, with a distance of 3–5 cm (1.2–2 inches).

4.Press the button and get the temperature.

Measure

7. Measure Surface Temperature

1.Press the measure button to turn on the machine.

Power on

2.Press the mode button to choose surface mode.

Choose mode

3.Point the probe at the center of surface, with a distance of 3–5 cm (1.2–2 inches).

4.Press the button and get the temperature.

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